

Anti-Human/Mouse/Rat ATG9A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04481

3th Edition

Description:Autophagy-related protein 9A is a protein that in humans is encoded by the ATG9A gene. Involved in autophagy and cytoplasm to vacuole transport (Cvt) vesicle formation. Plays a key role in the organization of the preautophagosomal structure/phagophore assembly site (PAS), the nucleating site for formation of the sequestering vesicle. Cycles between a juxta-nuclear trans-Golgi network compartment and late endosomes. Nutrient starvation induces accumulation on autophagosomes.

Antigen:Synthetic Peptide of Atg9a

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: