

**Anti-Human/Mouse/Rat COL4A1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04490**

---

3th Edition

**Description:** This gene encodes the major type IV alpha collagen chain of basement membranes. Like the other members of the type IV collagen gene family, this gene is organized in a head-to-head conformation with another type IV collagen gene so that each gene pair shares a common promoter.

**Antigen:** Synthesized peptide derived from Collagen IV

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Highly expressed in placenta.