

**Anti-Human FAM3D Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04495**

---

3th Edition

**Description:**FAM3D (Family With Sequence Similarity 3 Member D) is a Protein Coding gene. Diseases associated with FAM3D include Narcolepsy. GO annotations related to this gene include cytokine activity. An important paralog of this gene is FAM3A.

**Antigen:**Synthesized peptide derived from FAM3D

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Abundantly expressed in placenta and weakly expressed in small intestine.