

Anti-Human/Rat PLGRKT Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04512

3th Edition

Description: PLGRKT (Plasminogen Receptor With A C-Terminal Lysine) is a Protein Coding gene.

Antigen: Synthesized peptide derived from Plasminogen receptor

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in peripheral blood cells and monocytes. Expressed in adrenal medulla.