

**Anti-Human/Mouse/Rat SMARCB1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04526**

---

3th Edition

**Description:** The protein encoded by this gene is part of a complex that relieves repressive chromatin structures, allowing the transcriptional machinery to access its targets more effectively. The encoded nuclear protein may also bind to and enhance the DNA joining activity of HIV-1 integrase. This gene has been found to be a tumor suppressor, and mutations in it have been associated with Human, Mousealignant rhabdoid tumors. Two transcript variants encoding different isoforms have been found for this gene.

**Antigen:** Synthetic peptide from human protein

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**