

Anti-Human/Mouse/Rat IMP3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04527

3th Edition

Description: This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain. IMP3 (IMP3, U3 Small Nucleolar Ribonucleoprotein) is a Protein Coding gene. Diseases associated with IMP3 include Enchondroma. Among its related pathways are MAPK-Erk Pathway and rRNA processing in the nucleus and cytosol. GO annotations related to this gene include poly(A) RNA binding and structural constituent of ribosome.

Antigen: Synthetic peptide from human protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: