

Anti-Human/Mouse/Rat TPT1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04533

3th Edition

Description:Involved in calcium binding and microtubule stabilization.

Antigen:Synthesized peptide derived from the Internal region of human HRF

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Found in several healthy and tumoral cells including erythrocytes, hepatocytes, macrophages, platelets, keratinocytes, erythroleukemia cells, gliomas, melanomas, hepatoblastomas, and lymphomas. It cannot be detected in kidney and renal cell carcinoma (RCC). Expressed in placenta and prostate.