

Anti-Human/Mouse/Rat MYOZ2 Polyclonal Antibody**Polyclonal Antibody****Cat.NO.: PA04547**

3th Edition

Description:MYOZ2 (Myozenin 2) is a Protein Coding gene. Diseases associated with MYOZ2 include Cardiomyopathy, Hypertrophic, 16 and Myoz2-Related Familial Hypertrophic Cardiomyopathy. GO annotations related to this gene include actin binding and telethonin binding. An important paralog of this gene is MYOZ1. The protein encoded by this gene belongs to a family of sarcomeric proteins that bind to calcineurin, a phosphatase involved in calcium-dependent signal transduction in diverse cell types. These family members tether calcineurin to alpha-actinin at the z-line of the sarcomere of cardiac and skeletal muscle cells, and thus they are important for calcineurin signaling. Mutations in this gene cause cardiomyopathy familial hypertrophic type 16, a hereditary heart disorder.

Antigen:Synthetic peptide from human protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed specifically in heart and skeletal muscle.