

**Anti-Human/Mouse/Rat FAM48A Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04548**

---

3th Edition

**Description:**SUPT20H (SPT20 Homolog, SAGA Complex Component) is a Protein Coding gene. Among its related pathways are Autophagy - animal and Chromatin organization. GO annotations related to this gene include transcription cofactor activity.

**Antigen:**Synthetic peptide from human protein

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Highly expressed in testis, moderately in brain and pituitary gland. Expressed in several fetal tissues, including lung, brain, thymus and kidney. Expression is down-regulated in malignant prostate tissues.