

Anti-Human HSFX1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04551

3th Edition

Description:HSFX1 (Heat Shock Transcription Factor Family, X-Linked 1) is a Protein Coding gene. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and sequence-specific DNA binding. An important paralog of this gene is HSFX2.

Antigen:Synthesized peptide derived from the N-terminal region of human HSFX1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Testis-specific.