

Anti-Human/Mouse/Rat/Monkey HSPA5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04562

3th Edition

Description: Probably plays a role in facilitating the assembly of multimeric protein complexes inside the endoplasmic reticulum. Involved in the correct folding of proteins and degradation of misfolded proteins via its interaction with DNAJC10, probably to facilitate the release of DNAJC10 from its substrate.

Antigen: Synthesized peptide derived from the C-terminal region of human HSP A5

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: