

**Anti-Human/Mouse/Rat Histone H2A.Z Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04577**

---

3th Edition

**Description:** Histone H2A.Z/H2A.F/Z (H2A/z) is a 128 amino acid protein encoded by the human gene H2AFZ. Eukaryotic histones are basic and water soluble nuclear proteins that form hetero-octameric nucleosome particles by wrapping 146 base pairs of DNA sequentially in a left-handed super-helical turn to form chromosomal fiber. Two molecules of each of the four core histones (H2A, H2B, H3 and H4) form the octamer, which is comprised of two H2A-H2B dimers and two H3-H4 dimers, creating two nearly symmetrical halves by tertiary structure. H2A.Z/H2A.F/Z is a variant Histone H2A which replaces conventional H2A in a subset of nucleosomes. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of posttranslational modifications of histones, also called histone code, and nucleosome remodeling. H2A.Z/H2A.F/Z may be involved in the formation of constitutive heterochromatin and may be required for chromosome segregation during cell division.

**Antigen:** Synthetic peptide from human protein

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**