

**Anti-Human/Mouse/Rat HuB/D Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04582**

---

3th Edition

**Description:** Binds RNA. Seems to recognize a GAAA motif. Can bind to its own 3'-UTR, the FOS 3'-UTR and the ID 3'-UTR.

**Antigen:** Synthesized peptide derived from the Internal region of human HuB/D

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Brain.