

Anti-Human IL32 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04610

3th Edition

Description: Cytokine that may play a role in innate and adaptive immune responses. It induces various cytokines such as TNFA/TNF-alpha and IL8. It activates typical cytokine signal pathways of NF-kappa-B and p38 MAPK.

Antigen: Synthesized peptide derived from the C-terminal region of human IL-32.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Selectively expressed in lymphocytes. Expression more prominent in immune cells than in non-immune cells.