

Anti-Human ITPKC Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04616

3th Edition

Description: This gene encodes a member of the inositol 1,4,5-trisphosphate [Ins(1,4,5)P(3)] 3-kinase family of enzymes that catalyze the phosphorylation of inositol 1,4,5-trisphosphate to 1,3,4,5-tetrakisphosphate. The encoded protein is localized to the nucleus and cytoplasm and has both nuclear import and nuclear export activity. Single nucleotide polymorphisms in this gene are associated with Kawasaki disease. ITPKC (Inositol-Trisphosphate 3-Kinase C) is a Protein Coding gene. Diseases associated with ITPKC include Kawasaki Disease and Lymph Node Disease. Among its related pathways are Inositol phosphate metabolism (KEGG) and Calcium signaling pathway. GO annotations related to this gene include calmodulin binding and inositol-1,4,5-trisphosphate 3-kinase activity. An important paralog of this gene is ITPKB.

Antigen: Synthesized peptide derived from the Internal region of human InsP 3-kinase C

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in pancreas, skeletal muscle, liver, placenta and weakly in kidney and brain.