

Anti-Human IVL Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04620

3th Edition

Description:Involucrin, a component of the keratinocyte crosslinked envelope, is found in the cytoplasm and crosslinked to membrane proteins by transglutaminase. This gene is mapped to 1q21, among calpactin I light chain, trichohyalin, profillaggrin, loricrin, and calcyclin.

Antigen:Synthesized peptide derived from the C-terminal region of human Involucrin

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Keratinocytes of epidermis and other stratified squamous epithelia.