

**Anti-Human IRAK3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04628**

---

3th Edition

**Description:**Inhibits dissociation of IRAK1 and IRAK4 from the Toll-like receptor signaling complex by either inhibiting the phosphorylation of IRAK1 and IRAK4 or stabilizing the receptor complex.

**Antigen:**Synthesized peptide derived from the C-terminal region of human IRAK-M

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed predominantly in peripheral blood lymphocytes.