

Anti-Human/Mouse MAPK8IP3 Polyclonal Antibody**Polyclonal Antibody****Cat.NO.: PA04648**

3th Edition

Description:The protein encoded by this gene shares similarity with the product of Drosophila syd gene, required for the functional interaction of kinesin I with axonal cargo. Studies of the similar gene in mouse suggested that this protein may interact with, and regulate the activity of numerous protein kinases of the JNK signaling pathway, and thus function as a scaffold protein in neuronal cells. The C. elegans counterpart of this gene is found to regulate synaptic vesicle transport possibly by integrating JNK signaling and kinesin-1 transport. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. MAPK8IP3 (Mitogen-Activated Protein Kinase 8 Interacting Protein 3) is a Protein Coding gene. Among its related pathways are MAPK signaling pathway and Arf6 trafficking events. GO annotations related to this gene include kinesin binding and receptor signaling complex scaffold activity. An important paralog of this gene is SPAG9.

Antigen:Synthesized peptide derived from the Internal region of human JIP-3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: