

Anti-Human/Mouse/Rat PPP1R14C Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04663

3th Edition

Description:PP14C is an inhibitor of PPP1CA. It has over 600-fold higher inhibitory activity when phosphorylated, creating a molecular switch for regulating the phosphorylation status of PPP1CA substrates and smooth muscle contraction. It is found in breast carcinoma.

Antigen:Synthesized peptide derived from the Internal region of human KEPI

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Detected in breast cancer.