

Anti-Human/Mouse/Rat ZNF18 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04682

3th Edition

Description:ZNF18 (Zinc Finger Protein 18) is a Protein Coding gene. Among its related pathways are Gene Expression. GO annotations related to this gene include nucleic acid binding and RNA polymerase II transcription factor activity, sequence-specific DNA binding. An important paralog of this gene is ZNF202.

Antigen:Synthesized peptide derived from the Internal region of human KOX11

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: