

Anti-Human PMEPA1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04685

3th Edition

Description:Involved in down-regulation of the androgen receptor (AR). Enhances ubiquitination and proteasome-mediated degradation of AR, probably by recruiting NEDD4.

Antigen:Synthetic peptide from human protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highest expression in prostate. Also expressed in ovary.