

Anti-Human/Mouse RPA3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04696

3th Edition

Description:RPA3 (Replication Protein A3) is a Protein Coding gene. Diseases associated with RPA3 include Epidemic Typhus and Ovary Serous Adenocarcinoma. Among its related pathways are DNA Double-Strand Break Repair and Telomere C-strand (Lagging Strand) Synthesis. GO annotations related to this gene include single-stranded DNA binding and damaged DNA binding.

Antigen:Synthetic peptide from human protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: