

Anti-Human/Mouse/Rat ARHGEF2 Polyclonal Antibody
Polyclonal Antibody

Cat.NO.: PA04738

3th Edition

Description: Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form complex with G proteins and stimulate rho-dependent signals. Alternatively spliced transcript variants encoding different isoforms have been identified.

Antigen: Synthesized peptide derived from human Lfc around the non-phosphorylation site of Ser885.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: