



Anti-Human/Mouse/Rat LIMK2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04750

3th Edition

Description: Displays serine/threonine-specific phosphorylation of myelin basic protein and histone (MBP) in vitro.

Antigen: Synthesized peptide derived from the Internal region of human LIMK-2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highest expression in the placenta; moderate level in liver, lung, kidney, and pancreas. LIMK2a is found to be more abundant than LIMK2b in liver, colon, stomach, and spleen, while in brain, kidney, and placenta LIMK2b is the dominant form. In adult lung, both LIMK2a and LIMK2b is nearly equally observed.