

Anti-Human/Mouse STK11 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04762

3th Edition

Description:Essential role in G1 cell cycle arrest. Phosphorylates and activates members of the AMPK-related subfamily of protein kinases. Tumor suppressor.

Antigen:Synthesized peptide derived from human LKB1 around the non-phosphorylation site of Ser334.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitously expressed. Strongest expression in testis and fetal liver.