

Anti-Human/Mouse LMTK3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04767

3th Edition

Description:LMTK3 (Lemur Tyrosine Kinase 3) is a Protein Coding gene. Among its related pathways are GPCR Pathway and Nanog in Mammalian ESC Pluripotency. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is AATK.

Antigen:Synthesized peptide derived from the C-terminal region of human LMTK3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: