

Anti-Human/Mouse/Rat LRRC41 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04780

3th Edition

Description:LRRC41 (Leucine Rich Repeat Containing 41) is a Protein Coding gene. Among its related pathways are Innate Immune System and Class I MHC mediated antigen processing and presentation. GO annotations related to this gene include protein homodimerization activity.

Antigen:Synthesized peptide derived from the Internal region of human LRRC41

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: