

**Anti-Human/Mouse/Rat MAP3K13 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04794**

---

3th Edition

**Description:**The protein encoded by this gene is a member of serine/threonine protein kinase family. This kinase contains a dual leucine-zipper motif, and has been shown to form dimers/oligomers through its leucine-zipper motif. This kinase can phosphorylate and activate MAPK8/JNK, MAP2K7/MKK7, which suggests a role in the JNK signaling pathway.

**Antigen:**Synthesized peptide derived from the Internal region of human LZK

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in the adult brain, liver, placenta and pancreas, with expression strongest in the pancreas.