

Anti-Human/Mouse/Rat/Monkey MARCH3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04828

3th Edition

Description: This gene encodes a member of the membrane-associated RING-CH (MARCH) family. The encoded protein is an E3 ubiquitin-protein ligase that may be involved in regulation of the endosomal transport pathway. MARCH3 (Membrane Associated Ring-CH-Type Finger 3) is a Protein Coding gene. GO annotations related to this gene include ligase activity. An important paralog of this gene is MARCH2.

Antigen: Synthesized peptide derived from the Internal region of human March3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: