

Anti-Human/Mouse/Rat/Monkey MARK1/2/3/4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04832

3th Edition

Description: May play a role in cytoskeletal stability.

Antigen: Synthesized peptide derived from human MARK1/2/3/4 around the non-phosphorylation site of Thr215.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: MARK1 is Highly expressed in heart, skeletal muscle, brain, fetal brain and fetal kidney. MARK2 is High levels of expression in heart, brain, skeletal muscle and pancreas, lower levels observed in lung, liver and kidney.1 Publication MARK3 and MARK4 are Ubiquitous.