



Anti-Human/Mouse/Rat/Monkey MMP7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04848

3th Edition

Description:Degrades casein, gelatins of types I, III, IV, and V, and fibronectin. Activates procollagenase.

Antigen:Synthesized peptide derived from the C-terminal region of human MMP-7

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: