

Anti-Human/Mouse/Rat MC5R Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04865

3th Edition

Description: This gene encodes a member of the seven-pass transmembrane G protein-coupled melanocortin receptor protein family that stimulate cAMP signal transduction. The encoded protein is a receptor for melanocyte-stimulating hormone and adrenocorticotrophic hormone and is suggested to play a role in sebum generation.

Antigen: Synthesized peptide derived from the C-terminal region of human MC5-R

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in the brain but not in the melanoma cells.