

Anti-Human/Mouse/Rat MED18 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04892

3th Edition

Description: MED18 is a component of the Mediator complex, which is a coactivator for DNA-binding factors that activate transcription via RNA polymerase II (Sato et al., 2003 [PubMed 12584197]).

Antigen: Synthesized peptide derived from the Internal region of human Med18

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: