



**Anti-Human/Mouse MAP3K3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04901**

---

3th Edition

**Description:**Component of a protein kinase signal transduction cascade. Mediates activation of the NF-kappa-B, AP1 and DDIT3 transcriptional regulators.

**Antigen:**Synthesized peptide derived from the Internal region of human MEK Kinase-3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**