

Anti-Human/Mouse/Rat BMP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04903

3th Edition

Description:BMP1 (Bone Morphogenetic Protein 1) is a Protein Coding gene. Diseases associated with BMP1 include Osteogenesis Imperfecta, Type Xiii and Bmp1-Related Osteogenesis Imperfecta. Among its related pathways are GPCR Pathway and PEDF Induced Signaling. GO annotations related to this gene include calcium ion binding and growth factor activity. An important paralog of this gene is TLL2.This gene encodes a protein that is capable of inducing formation of cartilage in vivo. Although other bone morphogenetic proteins are members of the TGF-beta superfamily, this gene encodes a protein that is not closely related to other known growth factors. This gene is expressed as alternatively spliced variants that share an N-terminal protease domain but differ in their C-terminal region.

Antigen:Synthetic peptide from human protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous.