

**Anti-Human/Mouse CLEC3B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04911**

---

3th Edition

**Description:**CLEC3B (C-Type Lectin Domain Family 3 Member B) is a Protein Coding gene. Among its related pathways are Response to elevated platelet cytosolic Ca<sup>2+</sup>. GO annotations related to this gene include calcium ion binding and heparin binding. An important paralog of this gene is CLEC3A.

**Antigen:**Synthetic peptide from human protein

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Found in plasma.