



**Anti-Human MLANA Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04931**

---

3th Edition

**Description:**MLANA (Melan-A) is a Protein Coding gene. Diseases associated with MLANA include Melanoma and Balloon Cell Malignant Melanoma.

**Antigen:**Synthesized peptide derived from the Internal region of human Melan-A

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expression is restricted to melanoma and melanocyte cell lines and retina.