

Anti-Human CD52 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04977

3th Edition

Description:CD52 (CD52 Molecule) is a Protein Coding gene. Diseases associated with CD52 include Relapsing-Remitting Multiple Sclerosis and Chronic Lymphocytic Leukemia. Among its related pathways are Post-translational modification- synthesis of GPI-anchored proteins and Metabolism of proteins.

Antigen:Synthetic peptide from human protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: