

Anti-Human/Mouse/Rat MKP-5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05003

3th Edition

Description:Involved in the inactivation of MAP kinases. Has a specificity for the MAPK11/MAPK12/MAPK13/MAPK14 subfamily.

Antigen:Synthesized peptide derived from the C-terminal region of human MKP-5

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Detected in brain.