

Anti-Human/Mouse CNPY3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05010

3th Edition

Description:CNPY3 (Canopy FGF Signaling Regulator 3) is a Protein Coding gene. Diseases associated with CNPY3 include Swayback. Among its related pathways are Innate Immune System and Activated TLR4 signalling. GO annotations related to this gene include receptor binding. An important paralog of this gene is CNPY4. This gene encodes a protein that binds members of the toll-like receptor protein family and functions as a chaperone to aid in folding and export of these proteins. Alternative splicing results in multiple transcript variants. Naturally occurring readthrough transcription occurs between this locus and the downstream GNMT (glycine N-methyltransferase) gene and is represented with GenelD:107080644.

Antigen:Synthetic peptide from human protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: