

Anti-Human MMP13 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05019

3th Edition

Description:Degrades collagen type I. Does not act on gelatin or casein. Could have a role in tumoral process.

Antigen:Synthesized peptide derived from the N-terminal region of human MMP-13

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Seems to be specific to breast carcinomas.