

Anti-Human NUTM1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05038

3th Edition

Description:NUTM1 (NUT Midline Carcinoma Family Member 1) is a Protein Coding gene. Diseases associated with NUTM1 include Nut Midline Carcinoma and Extragonadal Nonseminomatous Germ Cell Tumor. Among its related pathways are Chromatin Regulation / Acetylation. An important paralog of this gene is NUTM2D.

Antigen:Synthesized peptide derived from human NUT

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Specifically expressed in testis.