

**Anti-Human/Mouse/Rat MRGPRF Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05055**

---

3th Edition

**Description:**MRGPRF (MAS Related GPR Family Member F) is a Protein Coding gene. GO annotations related to this gene include G-protein coupled receptor activity. An important paralog of this gene is MRGPRD.

**Antigen:**Synthesized peptide derived from the C-terminal region of human MRGF

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**