

Anti-Human MRGPRX4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05062

3th Edition

Description:MRGPRX4 (MAS Related GPR Family Member X4) is a Protein Coding gene. GO annotations related to this gene include G-protein coupled receptor activity. An important paralog of this gene is MRGPRX3.

Antigen:Synthesized peptide derived from the C-terminal region of human MRGX4

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Uniquely localized in a subset of small dorsal root and trigeminal sensory neurons.