

Anti-Human/Monkey ABCC10 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05068

3th Edition

Description:Involved in the transport of chloride ions. May regulate bicarbonate secretion and salvage in epithelial cells by regulating the SLC4A7 transporter.

Antigen:Synthesized peptide derived from the Internal region of human MRP7

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Isoform 1 is specifically expressed in spleen. Isoform 2 is more widely expressed.