

Anti-Human/Mouse/Rat MRPL14 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05080

3th Edition

Description: This nuclear gene encodes a protein component of the 39S subunit of the mitochondrial ribosome. A pseudogene of this gene is found on chromosome 17. MRPL14 (Mitochondrial Ribosomal Protein L14) is a Protein Coding gene. Among its related pathways are Mitochondrial translation and Viral mRNA Translation. GO annotations related to this gene include poly(A) RNA binding and structural constituent of ribosome.

Antigen: Synthesized peptide derived from the C-terminal region of human MRP-L14

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: