

Anti-Human OASL Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05116

3th Edition

Description: OASL encodes an interferon-induced protein that is a member of the 2-5A synthetase family. 2-prime,5-prime oligoadenylates (2-5As) bind to and activate RNase L, resulting in general RNA degradation and consequent inhibition of protein synthesis. 2-5As are produced by a well-conserved family of interferon-induced enzymes, the 2-5A synthetases or OASs

Antigen: Synthesized peptide derived from human OASL

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in most tissues, with the highest levels in primary blood Leukocytes and other hematopoietic system tissues, colon, stomach and to some extent in testis.