

Anti-Human/Mouse/Rat Phospho-EEF2 (Thr56) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05121

3th Edition

Description: This gene encodes a member of the GTP-binding translation elongation factor family. This protein is an essential factor for protein synthesis. It promotes the GTP-dependent translocation of the nascent protein chain from the A-site to the P-site of the ribosome. This protein is completely inactivated by EF-2 kinase phosphorylation.

Antigen: Synthesized phospho derived from human eEF2 (Phospho-Thr56)

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: