

Anti-Human/Mouse/Rat MRRF Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05123

3th Edition

Description:Responsible for the release of ribosomes from messenger RNA at the termination of protein biosynthesis. May increase the efficiency of translation by recycling ribosomes from one round of translation to another.

Antigen:Synthesized peptide derived from the C-terminal region of human MRRF

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: