

Anti-Human/Mouse/Rat APLF Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05165

3th Edition

Description: C2ORF13 is a component of the cellular response to chromosomal DNA single- and double-strand breaks. Nuclease involved in single-strand and double-strand DNA break repair. Recruited to sites of DNA damage through interaction with poly(ADP-ribose), a polymeric post-translational modification synthesized transiently at sites of chromosomal damage to accelerate DNA strand break repair reactions. Displays apurinic-aprimidinic (AP) endonuclease and 3'-5' exonuclease activities in vitro.

Antigen: Synthesized peptide derived from human APLF Polyclonal

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: